

# Ibex Transeat



The Ibex Transeat patient transport chair uses gravity and friction on an endless traction belt to transport the incapacitated up and down stairs, whilst its narrow track provides excellent manoeuvrability. It is easy to carry and simple to use, reducing the risk of back injury to ambulance and paramedic teams.

- Telescopic extendable handle for correct operator posture
- Stainless steel frame construction
- 159kg (25 stone) load capacity
- Compact storage with wall/floor mounted brackets
- Optional extras include dust cover and oxygen cylinder holder
- Adjustable patient straps
- Adjustable brake
- Quick and easy to open
- Detachable covers
- Floor & Wall Mounting Brackets

