
Rough Terrain Trolley











ParAid's Rough Terrain Trolley is a rugged stainless steel stretcher and patient transport trolley with wide pneumatic tyres for use in remote areas or over rough terrain, such as from helicopter landing pads to hospitals over gravel, grass and uneven ground. With hydraulically operated height adjustment of 800mm-1100mm, patient handling crews do not have to lift or carry the patient over distances, thus avoiding potential back injury and other musculoskeletal injuries. The ParAid Rough Terrain Trolley is designed using all-weather durable construction and materials for dependable outdoor use 365 days a year.

- Patient Load Capacity 200kg
- Flat, rubber covered stretcher platform 2000mm x 600mm
- Height adjustable hydraulic platform from 800mm to 1100mm
- 3 patient restraining belts to hold patient on stretcher/spine board
- Towing handle at each end

- 4 x 250mm dia pneumatic-tyred castors with brakes & direction locks
- Fold-out D/E size oxygen cylinder holder
- Telescopic IV pole with four hooks
- Foot brake and hand brake systems
- Fold-out monitor tray

